



COLD WATER PASTE GLAZE CONCENTRATE WM110

Review date 06/07/07

Macphie Cold Water Paste Glaze Concentrate is a concentrated paste which, when re-constituted with cold water, gives a rich, long-lasting, non-stick golden gloss when applied to all savoury pastry lines, bun goods, Danish pastries, etc.

RECIPE	kg	lb	oz	MAKE UP INSTRUCTIONS
Macphie Cold Water Paste Glaze Concentrate	1.000	2	4	<ul style="list-style-type: none">• Add the Macphie Cold Water Paste Glaze Concentrate to cold water and mix thoroughly until lump free• Allow the made up, liquid glaze to stand for 10 minutes and then re-mix briefly before use• The made up glaze should be applied by either brushing or spraying before baking• When using a spinning disc glazing machine Macphie Cold Water Paste Glaze Concentrate can be used at dilutions from 4-6 parts water to 1 part concentrate depending on the type of glaze required
Cold Water	5.000	11	4	
Total	6.000	13	8	

Please note:

A little separation in the pail may occur. If so, stir the mix in the pail before use if not using complete contents.

Please note:

Every care has been taken to ensure that the information detailed on this sheet is correct at the time of issue. However, your recipe, other ingredients and processing are outside our control. Therefore, we advise that you carry out prudent trials to check the product makes goods of a satisfactory quality and complies with current legislation.

Macphie of Glenberrie Ltd.,
Glenberrie, Stonehaven, Scotland, AB39 3YG. Tel +44 (0)1569 740641 Fax +44 (0)1569 740677
cservice@macphie.com www.macphie.com